



ENGLISH



Respect for Human Life and Pursuit of Beauty

ILWOO ELECTRON CO., LTD. is one of the leading home beauty and medical device(IPL) manufacturer in Korea.

We have been developing and manufacturing various kinds of products such as IPL device, Skin care device, and Hair care devices to provide quality products to our customers.

We also provide the differentiated OEM/ODM service with our ONE STOP SERVICE from the planning of the product, design, molding, injection and to the final mass production.

We will do our best to make the leap to the top global enterprise in the home beauty and medical device field.

CEO PARK, SOOJONG / ILWOO ELECTRON CO., LTD.

ABOUT US



COMPANY HISTORY

- 2024**
- Development and production of high-frequency massagers (ODM)
 - Development of medium frequency massager (ODM)
-

- 2020~2023**
- Renewal of ISO 13485:2016 Certificate
 - Renewal of GMP Conformity Certification
 - Patent Application No. 10-2114143 for Hand-piece Equipment for Intense Pulse Light (PL-3000)
 - Application for Skin Beauty Device Patent No. 10-2123667 with improved thermal conductivity and heat dissipation of Reevaa-S
 - Application for Patent No. 10-2123668 for Reevaa-S Composite Beauty Device
 - Development of Mass-Production Model for Skin Adhesive Hardening Machine
 - Development of model of PL-1500 (IPL)
 - Development and patent application of portable bimaxillary corrector (V-reevaa)
 - Development of Intense Pulse Light (IPL-1500)
 - Development and Patent of a Replaceable Scalp Massager (reevaa-H)

- 2018~2019**
- ISO 14001:2015 Certificate Renewal
 - Development of Intense Pulse Light (IPL-3000)
 - PL-3001 CVC Authentication
 - Change name of Promax Co., Ltd. to Ilwoo Electron Co., Ltd
 - Development of Fiona Beauty Device
 - Obtained CE, KC for SKIN SCRUBBER /FCC for reevaa-S, Sonic-H
 - Participate in a number of overseas exhibitions (e.g. Hong Kong COSMOPROF)
-

- 2007~2017**
- Obtained ISO 22765, ISO 14001 certification
 - Development of Intense Pulse Light and lamp cartridges
 - Development of Hair Iron, reevaa-H, reevaa-S,
 - Acquisition of KGMP Conformity Certificate, Obtained KC certificate, ISO 9001, ISO 14401, ISO 13485
 - Obtained Medical Device Manufacturing License by KFDA
 - Acquisition of KFDA medical device manufacturing license (PM-2640 SI/PM-2641 SI)
-

- 2006**
- Establishment of Promax Co., Ltd. (subsidiary company)
 - Development of Personal ultrasonic stimulator PM-2640 SI/PM-2641 SI (patent application)
-

- 1989**
- Establishment of ILWOO Precision Co., Ltd. (parent company)

Total Solution

From planning to design, production, and delivery, all the work is ONE-STOP PROCESS!

When the project begins, we allocate teams, share planned schedules, and share ideas that fit the product based on long-term know-how.

All departments participate in the project and proceed with an in-house consistent production system from planning/design/3D MOCK-Up-mold/injection/assembly production/quality inspection/delivery.

-
1. Production based on standardized work instructions.
 2. Quality improvement through inspection of first, second, and last items.
 3. Minimize risk by databasing defect occurrence information.
 4. The mold production period is 30 days from design to presentation.
-


INTEGRATED PRODUCTION SYSTEM



REVIEWING THE DESIGN PROPOSAL



INNER STRUCTURE AND EQUIPMENT DESIGN




PROTOTYPE CRITICISM AND SUPPLEMENTATION

OEM / ODM

20 years of accumulated technology and know-how!
IPL | ultrasonic | high frequency | medium frequency | galvanic
Development of medical and cosmetic devices



TECHNICAL MEETING AND PLAN



CORE AND SLIDE ASSEMBLY PRODUCT VISUAL EXAMINATION




SELECTION OF GOOD QUALITY PRODUCTS AND INSPECTION

CERTIFICATE

Through an insightful understanding of the product development process from planning and production of products, we carry out client-centered development and projects from woodworking, molds, and mass production.

We will do our best to become a reliable development partner who develops and supplies competitive products to clients by providing all the experiences and solutions ILWOO ELECTRON has.



38 Certificates and permits



Skin care & Hair removal

ANYMI IPL - 3000



Description

As a device that can take the strong light of a xenon lamp deep into the skin, make and treat various skin diseases.

It can be replaced with several treatments, such as hair removal, skin rejuvenation, and acne treatment at home.

The lifespan of the cartridge is a maximum 1 million times (based on multi-shot level 1).

A multi-shot function is added to this product for sensitive skin. This IPL product is certified by the Korean Ministry of Food and Drug Safety Association.

Cartridge replacement enables hair removal, whitening, and acne care.

How to use



Therapy

hair removal, skin rejuvenation, acne treatment.

Feature



Auto Shot



Multi Shot



Skin tone detection



Skin detection

Specification

Model name	ANYMI IPL - 3000
Voltage	AC 100 - 240V
Frequency	50/60hz
Weight	1.5kg
Product size	220 x 240 x 168 (mm)
Country of origin	Republic of Korea

Skin care & Hair removal

ANYMI IPL - 2000



Description

As a device that can take the strong light of a xenon lamp deep into the skin and make and treat various skin diseases.


It can be replaced with several treatments, such as hair removal, skin rejuvenation, and acne treatment at home.

The lifespan of the cartridge is a maximum 150,000 times (based on level 1).

This IPL product is certified by the Korean Ministry of Food and Drug Safety Association.

Cartridge replacement enables hair removal, whitening, and acne care.

How to use



Therapy

hair removal, skin rejuvenation, acne treatment.

Feature



Auto Shot



Skin detection

Specification

Model name	ANYMI IPL - 2000
Voltage	AC 100 - 240V
Frequency	50/60hz
Weight	1.5kg
Product size	236 x 259 x 167 (mm)
Country of origin	Republic of Korea

Skin care & Hair removal

ANYMI IPL-1500



Description

As a device that can take the strong light of a xenon lamp deep into the skin, make and treat various skin diseases,

It can be replaced with several treatments, such as hair removal, skin rejuvenation, and acne treatment at home.

The lifespan of the cartridge is a maximum 100,000 times (based on level 1).

This IPL product is certified by the Korean Ministry of Food and Drug Safety Association.

It is available 1,000,000 times (based on multi-level 1), and this product has added a multi-pulse function that surveys 4 times in a row for a short time.

Cartridge replacement enables hair removal, whitening, and acne care.

How to use



Therapy

hair removal, skin rejuvenation, acne treatment.

Feature



Auto Shot



Skin detection



4S Multi Shot

사양

Model name	ANYMI IPL - 1500 (ANYMI PLUS)
Voltage	AC 100 - 240V
Frequency	50/60Hz
Weight	1.5kg
Product size	196 x 139 x 180 (mm)
Country of origin	Republic of Korea

Catridge

Skin detection Sensor

The skin detection sensor is a safety device to protect the user from pulsed light and prevent the skin from the IPL shot if the skin is not detected by the skin detection sensor.

APPLIED PRODUCTS

ANIMY : IPL - 3000, IPL - 2000, IPL - 1000

Skin Tone Sensor

The skin tone sensor instantly adjusts the intensity of the light according to the user's skin tone for skin safety.

APPLIED PRODUCTS


ANIMY IPL - 3000



ANIMY IPL - 3000



ANIMY IPL - 2000



ANIMY IPL - 1500



Classification (based on ANIMY IPL 3000)

	640nm H-type cartridge (for Hair Removal)	580nm -type cartridge (for Skin Rejuvenation)	380nm A-type cartridge (for Acne Treatment)
Pulse wave type	square wave	square wave	square wave
Pulse Width	7ms ± 10%	7ms ± 10%	7ms ± 10%
Multi shot	1 pulse	1 pulse	1 pulse
Output stability	±20%	±20%	±20%
Shot Area	32mm×22mm±10%	32mm×22mm±10%	32mm×22mm±10%
Output Energy Adjustment Unit	0.10 ~ 1.25 J/cm ²	0.10 ~ 1.25 J/cm ²	0.10 ~ 1.25 J/cm ²
Output Wavelength	640 ~ 900 nm ± 10%	550~900nm ± 10%	380 ~ 700 nm ± 10%

Hot & Cool For Skin care
reevaa S



Description

A massager capable of galvanic and ultrasonic heating and cooling functions as a single device that absorbs the active ingredients of cosmetics deep into the skin layer or, conversely. Discharges waste from the skin layer by using the principle of pushing and pulling positive and negative ions of galvanic ions.

How to use



Therapy

Galvanic / ultrasonic / thermal / cooling functions are all in one device. Clean mode, Nutrition mode, and Boost mode combine the functions required for each mode smartly. 3 MHz ultrasound, which vibrates 3 million times per second, gives hypoallergenic massage to sensitive facial skin. The cooling function calms the skin and helps with elasticity. Adjustable output, ion strength, and temperature. Possession of patented technology. Automatic off function after 12 minutes for skin stabilization.

Feature

-  Galvanic ion
-  Ultra sonic
-  Cleansing Massage
-  Nutrition Massage
-  Lifting Massage
-  Hot & Cold Massage

Specification


Model name	reevaa S
Voltage	3.6V 2600mA Lithium ion battery
Weight	207g
Product size	42 x 65 x 179 (mm)
Country of origin	Republic of Korea

Feature by Mode

1. CLEAN MODE

- The heat at approximately 42°C opens skin pores and improves blood circulation.
- ⊕ion pulls waste products to clean skin and exfoliate skin.

Differences may occur depending on the usage environment and settings.




2. NUTRITION MODE

- The heat at approximately 42°C opens skin pores and improves blood circulation.
- Microvibration 3 million times per second with 3 mHz ultrasonic waves provides hypoallergenic massage and stimulates cells in the skin to improve elasticity.
- ⊖ions penetrate the active ingredients of cosmetics deeply.



3. BOOSTING MODE

- The cooling function at approximately 12°C closes pores and soothes the skin.
- ⊕⊖ ions come out at the same time to give skin elasticity.



Ultrasonic, far infrared, ion for body therapy

reevaa **FIT**



Description

REEVAA FIT with warm ultrasonic waves. REEVAA FIT ultrasonic thermal massager absorbs 500 kHz ultrasonic waves and galvanics to the maximum of the main ingredients of the body massage cream or gel. Thermal mode, which emits far-infrared rays, helps improve circulation. This stimulates the cellulite that clumps up in the skin, releasing it smoothly and creating a firm body line. It also helps muscles that need a thermal compress.

How to use



Clinically proven effectiveness

8 weeks
after use

Cellulite reduction **29.44%**
Elasticity improvement **15.30%**
Blood circulation improvement **11.34%**
Relief of swelling **11.20%**

(Test agency: SKINMED Co., Ltd., 2020.07.12~2023.09.08, human body application test report on the reduction of REEVAA FIT cellulite, human body application test report on the improvement of REEVAA FIT elasticity, and human body application test report on the parallel use of REEVAA FIT PPC elasticity cream) ※ There may be individual differences.

Feature



Ultra
sonic



Galvanic
ion



Far
infrared



Vibration

Specification

Model name	reevaa fit
Power	DC 24V 2.5A
Weight	259g
Product size	120 x 11 x 101 (mm)
Country of origin	Republic of Korea

How to take care of REEVAAFIT

1. Elastic and healthy body line care

Use Body Care Mode + Thermal Mode Step 1

- LOOSE BELLY FAT
- HARD CELLULITE
- BUMPY FLESH
- LOOSE FLESH



500 kHz powerful ultrasonic waves penetrate deep into the skin to stimulate cells, and thermal far-infrared rays help blood circulation to stimulate muscles and cellulite clumped in the skin, softening and creating a solid body line.



Differences may occur depending on the usage environment and settings.


2. Uncomfortable muscle care

Use Body Care Mode + Thermal Mode Step 2

- CLUMPED MUSCLES
- A STIFF SHOULDER
- WHEN YOUR WRISTS, ANKLES, AND KNEES ARE UNCOMFORTABLE
- WHEN YOUR LEGS ARE SWOLLEN



500 kHz powerful ultrasonic waves release hardened muscles and thermal far-infrared rays help blood circulation and improve muscle relaxation, muscle cramps, and joint flexibility to relieve muscle stress.



3. Body circulation care

Thermal Mode Step 2

- Where blood circulation care is needed
- When my stomach hurts because of menstrual cramps
- When your hands and feet are cold
- When the body temperature is low
- When you feel uncomfortable with varicose veins.



The heat and ceramic far-infrared rays are deep inside the skin. It transfers heat to strengthen blood flow and help the body circulate.



Scalp massage reevaa H⁺



Description

Two brushes that help with hair and scalp care can be used conveniently depending on the application. Ball bearing brush is a scalp care beauty device that maximizes scalp care effects by applying nutrients evenly to help absorb useful ingredients well.

How to use



Therapy

If you use the Hair Growing Redients product after inserting it into the ball bearing brush, it applies hair growing redients to the scalp hygienically and conveniently, and maximizes the effect of LED infrared rays, galvanic ions negative (-) ions, and vibrations on the scalp.

Feature



Galvanic ion



infrared



Vibration Massage



absorption

Specification

Model name	reevaa H+
Voltage	3.7V 250mA Lithium ion battery
Weight	65g
Product size	211 x 75 x 20 (mm)
Country of origin	Republic of Korea


Features and Strengths of REVAAH+

- Ball bearing application brush
- With the introduction of galvanic ions, drugs are injected deep into the scalp
- LED infrared wavelengths transmit light to subcutaneous tissue to activate cells
- Up to 10,000 microvibration scalp massages per minute
- Hygienic management with removable replaceable brushes
- Easy charging with C-Type charging jack



Features of the Drug Application Brush

- Apply evenly throughout the scalp
- Effectively absorbs deep into the scalp barrier with galvanic \ominus ions and vibrations.



Features of LED brushes

- Transmit light and heat to the upper layer of the dermis.
- Stimulating scalp cells and strengthening scalp barriers.



Scalp massage and cleansing

ANYMI SONIC H (HAIR SONIC)



Description

It is a home hair care massager that helps useful hair product ingredients absorbed by using ultrasonic and galvanic ions.

UltraSonic stimulates the scalp to mitigate, the amount (+) of galvanic ions helps clear the scalp, and the useful ingredients of negative (-) ions and woolen agents are absorbed in the scalp to maximize the effect.

How to use



Therapy

Ultrasound and micro-vibration make your scalp healthier.



Galvanic ion



Ultra sonic



Cleansing Massage



Nutrition Massage



Vibration Massage

Specification

Model name	ANYMI SONIC H (HAIR SONIC)
Voltage	3.7V 1000mA Lithium ion battery
Weight	185g
Product size	185 x 81 x 38 (mm)
Country of origin	Republic of Korea

Other Developed Products



PRO IRON



MINI IRON



TWEEZER



SCRUBBER



LED ION

Other Developed Products



SONIC F



HAIR ION



BB UP




4 CH



V REEVAA



"Company Realizes the Impression of Customers"



ILWOOELECTRON CO.,LTD.
www.ilwooelectron.com

[Factory in Bucheon] 21, Saneop-ro 7beon-gil, Ojeong-gu, Bucheon-si, Gyeonggi-do, Republic of Korea
[Factory in Incheon] 72, Seounsandan-ro 2-gil, Gyeyang-gu, Incheon, Republic of Korea

TEL +82.2.839.2011 | FAX +82.2.839.4782 | E - mail ilwooelectron@ilwooelectron.com